



PSMB9 Antibody, Biotin conjugated

Product Code	CSB-PA018887LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28065
Immunogen	Recombinant Human Proteasome subunit beta type-9 protein (21-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit beta type-9, Low molecular mass protein 2, Macropain chain 7, Multicatalytic endopeptidase complex chain 7, Proteasome chain 7, Proteasome subunit beta-1i, Really interesting new gene 12 protein, PSMB9, LMP2, PSMB6i, RING12
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	PSMB9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.