




VAMP4 Antibody

Product Code	CSB-PA025783LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75379
Immunogen	Recombinant Human Vesicle-associated membrane protein 4 protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle-associated membrane protein 4, VAMP-4, VAMP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	VAMP4

Image



Western Blot

Positive WB detected in: K562 whole cell lysate, RAW264.7 whole cell lysate

All lanes: VAMP4 antibody at 4.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 17 kDa

Observed band size: 17 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.