



CFHR1 Antibody, Biotin conjugated

Product Code	CSB-PA13809D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03591
Immunogen	Recombinant Human Complement factor H-related protein 1 protein (132-221AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement factor H-related protein 1, FHR-1, H factor-like protein 1, H-factor-like 1, H36, CFHR1, CFHL, CFHL1, CFHL1P, CFHR1P, FHR1, HFL1, HFL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular; Immunology
Target Names	CFHR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.