



KANK2 Antibody, Biotin conjugated

Product Code	CSB-PA720638LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63ZY3
Immunogen	Recombinant Human KN motif and ankyrin repeat domain-containing protein 2 protein (410-585AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	KN motif and ankyrin repeat domain-containing protein 2, Ankyrin repeat domain-containing protein 25, Matrix-remodeling-associated protein 3, SRC-1-interacting protein, SIP, SRC-interacting protein, SRC1-interacting protein, KANK2, ANKRD25, KIAA1518, MXRA3, SIP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KANK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.