



# STRA8 Antibody, Biotin conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA767212LD01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                            |
| Uniprot No.         | Q7Z7C7   |
| Immunogen           | Recombinant Human Stimulated by retinoic acid gene 8 protein homolog protein (115-241AA) |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Biotin   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4         |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Stimulated by retinoic acid gene 8 protein homolog, STRA8                                |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Developmental biology; Stem cells  |
| Target Names        | STRA8  |