



PRPF39 Antibody, HRP conjugated

Product Code	CSB-PA769774LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UA1
Immunogen	Recombinant Human Pre-mRNA-processing factor 39 protein (41-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA-processing factor 39, PRP39 homolog, PRPF39
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRPF39
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.