



PHLDB1 Antibody, Biotin conjugated

Product Code	CSB-PA773033LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UU1
Immunogen	Recombinant Human Pleckstrin homology-like domain family B member 1 protein (233-423AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology-like domain family B member 1, Protein LL5-alpha, PHLDB1, KIAA0638, LL5A, DLNB07
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	PHLDB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.