



NSL1 Antibody, HRP conjugated

Product Code	CSB-PA839345LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IY1
Immunogen	Recombinant Human Kinetochore-associated protein NSL1 homolog protein (225-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinetochore-associated protein NSL1 homolog, NSL1, C1orf48, DC31, DC8, MIS14
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NSL1