



BARD1 Antibody, Biotin conjugated

Product Code	CSB-PA860351LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99728
Immunogen	Recombinant Human BRCA1-associated RING domain protein 1 protein (428-777AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BRCA1-associated RING domain protein 1, RING-type E3 ubiquitin transferase BARD1, , BARD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling; Cancer; Cell biology; Metabolism
Target Names	BARD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.