



EFHD1 Antibody, FITC conjugated

Product Code	CSB-PA880107LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUP0
Immunogen	Recombinant Human EF-hand domain-containing protein D1 protein (15-72AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EF-hand domain-containing protein D1, EF-hand domain-containing protein 1, Swiprosin-2, EFHD1, SWS2, PP3051
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	EFHD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.