



LIME1 Antibody, Biotin conjugated

Product Code	CSB-PA884463LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H400
Immunogen	Recombinant Human Lck-interacting transmembrane adapter 1 protein (170-295AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lck-interacting transmembrane adapter 1, Lck-interacting membrane protein, Lck-interacting molecule, LIME1, LIME, LP8067
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LIME1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.