



SPINK5 Antibody, HRP conjugated

Product Code	CSB-PA889077LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ38
Immunogen	Recombinant Human Serine protease inhibitor Kazal-type 5 protein (722-847AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor Kazal-type 5, Lympho-epithelial Kazal-type-related inhibitor, LEKTI, , Hemofiltrate peptide HF6478, Hemofiltrate peptide HF7665, SPINK5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	SPINK5