



SPINK5 Antibody, Biotin conjugated

Product Code	CSB-PA889077LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ38
Immunogen	Recombinant Human Serine protease inhibitor Kazal-type 5 protein (722-847AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor Kazal-type 5, Lympho-epithelial Kazal-type-related inhibitor, LEKTI, , Hemofiltrate peptide HF6478, Hemofiltrate peptide HF7665, SPINK5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	SPINK5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.