



# PLAGL1 Antibody, Biotin conjugated

Product Code	CSB-PA892351LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM63
Immunogen	Recombinant Human Zinc finger protein PLAGL1 protein (237-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein PLAGL1, Lost on transformation 1, LOT-1, Pleiomorphic adenoma-like protein 1, Tumor suppressor ZAC, PLAGL1, LOT1, ZAC
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling; Neuroscience; Stem cells
Target Names	PLAGL1