



LSM11 Antibody, FITC conjugated

Product Code	CSB-PA013201LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83369
Immunogen	Recombinant Human U7 snRNA-associated Sm-like protein LSm11 protein (1-141AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U7 snRNA-associated Sm-like protein LSm11, LSM11
Immunogen Species	Homo sapiens (Human)
Target Names	LSM11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.