



PDE8B Antibody, Biotin conjugated

Product Code	CSB-PA529743LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95263
Immunogen	Recombinant Human High affinity cAMP-specific and IBMX-insensitive 3',5'-cyclic phosphodiesterase 8Bprotein (18-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity cAMP-specific and IBMX-insensitive 3',5'-cyclic phosphodiesterase 8B (HsPDE8B) (EC 3.1.4.53) (Cell proliferation-inducing gene 22 protein), PDE8B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	PDE8B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.