




TTC39B Antibody

Product Code	CSB-PA713348LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VTQ0
Immunogen	Recombinant Human Tetratricopeptide repeat protein 39Bprotein (10-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 39B (TPR repeat protein 39B), TTC39B, C9orf52
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	TTC39B

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, MCF-7 whole cell lysate

All lanes: TTC39B antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 77, 67, 76, 60, 70, 23 kDa

Observed band size: 77 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.