



TTC39B Antibody, HRP conjugated

Product Code	CSB-PA713348LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VTQ0
Immunogen	Recombinant Human Tetratricopeptide repeat protein 39Bprotein (10-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 39B (TPR repeat protein 39B), TTC39B, C9orf52
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	TTC39B