



TTC39B Antibody, Biotin conjugated

Product Code	CSB-PA713348LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VTQ0
Immunogen	Recombinant Human Tetratricopeptide repeat protein 39Bprotein (10-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 39B (TPR repeat protein 39B), TTC39B, C9orf52
Immunogen Species	Homo sapiens (Human)
Research Area	Cell biology
Target Names	TTC39B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.