



DPH5 Antibody, Biotin conjugated

Product Code	CSB-PA862001LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2P9
Immunogen	Recombinant Human Diphthine methyl ester synthase protein (138-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphthine methyl ester synthase (EC 2.1.1.314) (Diphthamide biosynthesis methyltransferase), DPH5
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer; Metabolism; Signal transduction
Target Names	DPH5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.