



GRHL1 Antibody, FITC conjugated

Product Code	CSB-PA868368LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZI5
Immunogen	Recombinant Human Grainyhead-like protein 1 homologprotein (1-95AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Grainyhead-like protein 1 homolog (Mammalian grainyhead) (NH32) (Transcription factor CP2-like 2) (Transcription factor LBP-32), GRHL1, LBP32 MGR TFPCP2L2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GRHL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.