



VPS51 Antibody, HRP conjugated

Product Code	CSB-PA883424LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UID3
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 51 homolog protein (595-673AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 51 homolog (Another new gene 2 protein) (Protein fat-free homolog), VPS51, ANG2 C11orf2 C11orf3 FFR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	VPS51