



VPS51 Antibody, Biotin conjugated

Product Code	CSB-PA883424LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UID3
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 51 homolog protein (595-673AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 51 homolog (Another new gene 2 protein) (Protein fat-free homolog), VPS51, ANG2 C11orf2 C11orf3 FFR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal transduction
Target Names	VPS51
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.