




ORF7 Antibody

Product Code	CSB-PA21639A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	V5N6H2
Immunogen	Recombinant Porcine reproductive and respiratory syndrome virus Nucleocapsid protein (1-123AA)
Raised In	Rabbit
Species Reactivity	Porcine reproductive and respiratory syndrome virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ORF7
Immunogen Species	Porcine reproductive and respiratory syndrome virus
Target Names	ORF7

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ORF7 antibody at 2.1µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.