



ORF7 Antibody, Biotin conjugated

Product Code	CSB-PA21639D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	V5N6H2
Immunogen	Recombinant Porcine reproductive and respiratory syndrome virus Nucleocapsid protein (1-123AA)
Raised In	Rabbit
Species Reactivity	Porcine reproductive and respiratory syndrome virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ORF7
Immunogen Species	Porcine reproductive and respiratory syndrome virus
Target Names	ORF7