



ORF7 Antibody, Biotin conjugated

Product Code	CSB-PA21639D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	V5N6H2
Immunogen	Recombinant Porcine reproductive and respiratory syndrome virus Nucleocapsid protein (1-123AA)
Raised In	Rabbit
Species Reactivity	Porcine reproductive and respiratory syndrome virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ORF7
Immunogen Species	Porcine reproductive and respiratory syndrome virus
Target Names	ORF7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.