




# cas3 Antibody

<b>Product Code</b>	CSB-PA24029A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q53VY2
<b>Immunogen</b>	Recombinant Thermus thermophilus CRISPR-associated endonuclease/helicase Cas3 protein (120-226AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Thermus thermophilus
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CRISPR-associated endonuclease/helicase Cas3 (EC 3.1.-.-) (EC 3.6.4.-) (CRISPR-associated ssDNA endonuclease/helicase Cas3), cas3
<b>Immunogen Species</b>	Thermus thermophilus
<b>Target Names</b>	cas3

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: cas3 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.