



cas3 Antibody, Biotin conjugated

Product Code	CSB-PA24029D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53VY2
Immunogen	Recombinant Thermus thermophilus CRISPR-associated endonuclease/helicase Cas3 protein (120-226AA)
Raised In	Rabbit
Species Reactivity	Thermus thermophilus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CRISPR-associated endonuclease/helicase Cas3 (EC 3.1.-.-) (EC 3.6.4.-) (CRISPR-associated ssDNA endonuclease/helicase Cas3), cas3
Immunogen Species	Thermus thermophilus
Target Names	cas3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.