




skib Antibody

Product Code	CSB-PA24359A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9YI02
Immunogen	Recombinant Zebrafish Nuclear oncoprotein protein (289-438AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	skib
Immunogen Species	Zebrafish
Target Names	skib

Image



Western Blot

Positive WB detected in: Zebrafish brain tissue
 All lanes: skib antibody at 5µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 80 kDa
 Observed band size: 80 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.