



SDR Antibody, FITC conjugated

Product Code	CSB-PA24399C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A0A0
Immunogen	Recombinant Zebrafish Short-chain alcohol dehydrogenase protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SDR
Immunogen Species	Zebrafish
Target Names	SDR