




# src Antibody

<b>Product Code</b>	CSB-PA24479A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q1JPZ3
<b>Immunogen</b>	Recombinant Zebrafish Proto-oncogene tyrosine-protein kinase Src protein (2-146AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Zebrafish
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Proto-oncogene tyrosine-protein kinase Src (EC 2.7.10.2) (Proto-oncogene c-Src) (pp60c-src) (p60-Src), src
<b>Immunogen Species</b>	Zebrafish
<b>Target Names</b>	src

## Image



### Western Blot

Positive WB detected in: Zebrafish tissue

All lanes: src antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 61 kDa

Observed band size: 61 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.