



pepP Antibody, HRP conjugated

Product Code	CSB-PA303407LB01MLW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P75313
Immunogen	Recombinant Mycoplasma pneumoniae Putative Xaa-Pro aminopeptidase protein (1-354AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative Xaa-Pro aminopeptidase (X-Pro aminopeptidase) (EC 3.4.11.9) (Aminoacylproline aminopeptidase) (Aminopeptidase P) (APP), pepP
Immunogen Species	Mycoplasma pneumoniae
Target Names	pepP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.