



MPN_083 Antibody, Biotin conjugated

Product Code	CSB-PA303444LD01Mlw
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P75610
Immunogen	Recombinant Mycoplasma pneumoniae Uncharacterized lipoprotein MPN_083 protein (22-242AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized lipoprotein MPN_083
Immunogen Species	Mycoplasma pneumoniae
Target Names	MPN_083
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.