




PRA1 Antibody

Product Code	CSB-PA308545YA01CZD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87020
Immunogen	Recombinant Candida albicans pH-regulated antigen PRA1 protein (16-299AA)
Raised In	Rabbit
Species Reactivity	Candida albicans
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	pH-regulated antigen PRA1 (58 kDa fibrinogen-binding mannoprotein), PRA1, FBP1
Immunogen Species	Candida albicans
Target Names	PRA1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: PRA1 antibody at 4.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.