



PRA1 Antibody, Biotin conjugated

Product Code	CSB-PA308545YD01CZD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87020
Immunogen	Recombinant Candida albicans pH-regulated antigen PRA1 protein (16-299AA)
Raised In	Rabbit
Species Reactivity	Candida albicans
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	pH-regulated antigen PRA1 (58 kDa fibrinogen-binding mannoprotein), PRA1, FBP1
Immunogen Species	Candida albicans
Target Names	PRA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.