




BHRF1 Antibody

Product Code	CSB-PA314488LA01EFC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03182
Immunogen	Recombinant Epstein-Barr virus Apoptosis regulator BHRF1 protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis regulator BHRF1 (Early antigen protein R) (EA-R) (Nuclear antigen)
Immunogen Species	Epstein-Barr virus
Target Names	BHRF1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: BHRF1 antibody at 4μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa