




oprF Antibody

Product Code	CSB-PA318417LA01EZ X
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13794
Immunogen	Recombinant Pseudomonas aeruginosa Outer membrane porin F protein (25-350AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane porin F, oprF
Immunogen Species	Pseudomonas aeruginosa
Target Names	oprF

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: oprF antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39.3 kDa

Observed band size: 44 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.