



oprF Antibody, HRP conjugated

Product Code	CSB-PA318417LB01EZ X
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13794
Immunogen	Recombinant Pseudomonas aeruginosa Outer membrane porin F protein (25-350AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane porin F, oprF
Immunogen Species	Pseudomonas aeruginosa
Target Names	oprF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.