




rsrIM Antibody

Product Code	CSB-PA318583LA01RLE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14751
Immunogen	Recombinant Rhodobacter sphaeroides Modification methylase RsrI protein (1-319AA)
Raised In	Rabbit
Species Reactivity	Rhodobacter sphaeroides
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Modification methylase RsrI (M.RsrI) (EC 2.1.1.72) (Adenine-specific methyltransferase RsrI), rsrIM
Immunogen Species	Rhodobacter sphaeroides
Target Names	rsrIM

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: rsrIM antibody at 4.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.