




BLLF3 Antibody

Product Code	CSB-PA319032LA01EFC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	K9US42
Immunogen	Recombinant Epstein-Barr virus Deoxyuridine 5'-triphosphate nucleotidohydrolase protein (1-278AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyuridine 5'-triphosphate nucleotidohydrolase, dUTP pyrophosphatase, BLLF3
Immunogen Species	Epstein-Barr virus
Target Names	BLLF3

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: BLLF3 antibody at 3.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.