



# BLLF3 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA319032LB01EFC
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	K9US42
<b>Immunogen</b>	Recombinant Epstein-Barr virus Deoxyuridine 5'-triphosphate nucleotidohydrolase protein (1-278AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Epstein-Barr virus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Deoxyuridine 5'-triphosphate nucleotidohydrolase, dUTP pyrophosphatase, BLLF3
<b>Immunogen Species</b>	Epstein-Barr virus
<b>Target Names</b>	BLLF3