



BLLF3 Antibody, Biotin conjugated

Product Code	CSB-PA319032LD01EFC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	K9US42
Immunogen	Recombinant Epstein-Barr virus Deoxyuridine 5'-triphosphate nucleotidohydrolase protein (1-278AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyuridine 5'-triphosphate nucleotidohydrolase, dUTP pyrophosphatase, BLLF3
Immunogen Species	Epstein-Barr virus
Target Names	BLLF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.