



# S Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA322803LC01BJK
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P15777
<b>Immunogen</b>	Recombinant Bovine coronavirus Spike glycoprotein protein (326-540AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bovine coronavirus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Spike glycoprotein (S glycoprotein) (E2) (Peplomer protein) [Cleaved into: Spike protein S1 (90B); Spike protein S2 (90A)], S
<b>Immunogen Species</b>	Bovine coronavirus
<b>Target Names</b>	S
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.