



S Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA322803LD01BJK |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P15777 |
| Immunogen | Recombinant Bovine coronavirus Spike glycoprotein protein (326-540AA) |
| Raised In | Rabbit |
| Species Reactivity | Bovine coronavirus |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Spike glycoprotein (S glycoprotein) (E2) (Peplomer protein) [Cleaved into: Spike protein S1 (90B); Spike protein S2 (90A)], S |
| Immunogen Species | Bovine coronavirus |
| Target Names | S |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |