



vvhA Antibody, Biotin conjugated

Product Code	CSB-PA323085LD01VFI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19247
Immunogen	Recombinant Vibrio vulnificus Cytolysin protein (21-471AA)
Raised In	Rabbit
Species Reactivity	Vibrio vulnificus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytolysin, vvhA
Immunogen Species	Vibrio vulnificus
Target Names	vvhA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.