




C Antibody

Product Code	CSB-PA323345LA01RAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15201
Immunogen	Recombinant Rabies virus Putative uncharacterized protein C protein (1-69AA)
Raised In	Rabbit
Species Reactivity	Rabies virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein C, C
Immunogen Species	Rabies virus
Target Names	C

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: C antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.