




fbpB Antibody

Product Code	CSB-PA324138LA01MKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21160
Immunogen	Recombinant Mycobacterium kansasii Diacylglycerol acyltransferase/mycolyltransferase Ag85B protein (41-325AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium kansasii
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol acyltransferase/mycolyltransferase Ag85B (DGAT) (EC 2.3.1.122) (EC 2.3.1.20) (30 kDa extracellular protein) (Acyl-CoA: diacylglycerol acyltransferase) (Antigen 85 complex B) (85B) (Ag85B) (Extracellular alpha-antigen) (Fibronectin-binding protein B) (Fbps B), fbpB
Immunogen Species	Mycobacterium kansasii
Target Names	fbpB

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: fbpB antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.