



# fbpB Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA324138LD01MKZ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P21160
<b>Immunogen</b>	Recombinant Mycobacterium kansasii Diacylglycerol acyltransferase/mycolyltransferase Ag85B protein (41-325AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mycobacterium kansasii
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Diacylglycerol acyltransferase/mycolyltransferase Ag85B (DGAT) (EC 2.3.1.122) (EC 2.3.1.20) (30 kDa extracellular protein) (Acyl-CoA: diacylglycerol acyltransferase) (Antigen 85 complex B) (85B) (Ag85B) (Extracellular alpha-antigen) (Fibronectin-binding protein B) (Fbps B), fbpB
<b>Immunogen Species</b>	Mycobacterium kansasii
<b>Target Names</b>	fbpB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.