



ompA Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA332793LB01HSZ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P45996 |
| Immunogen | Recombinant Haemophilus influenzae Outer membrane protein P5 protein (22-359AA) |
| Raised In | Rabbit |
| Species Reactivity | Haemophilus influenzae |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Outer membrane protein P5 (OMP P5) (Fimbrin), ompA, ompP5 |
| Immunogen Species | Haemophilus influenzae |
| Target Names | ompA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |