



ompA Antibody, FITC conjugated

Product Code	CSB-PA332793LC01HSZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45996
Immunogen	Recombinant Haemophilus influenzae Outer membrane protein P5 protein (22-359AA)
Raised In	Rabbit
Species Reactivity	Haemophilus influenzae
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein P5 (OMP P5) (Fimbrin), ompA, ompP5
Immunogen Species	Haemophilus influenzae
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.